

## LEGENDA








### N.B. - Codici colori - alogene/al tungsteno:









Lampade alogene 2,8V hanno involucro dorato con base blu.  
Lampade alogene 3,5V hanno involucro dorato con base rossa.  
Lampade al tungsteno 2,8V hanno involucro argentato con base blu.

### N.B. - Colour Codes - Halogen/Tungsten:

Halogen 2.8V bulbs have gold cases with blue base.  
Halogen 3.5V bulbs have gold cases with red base.  
Tungsten 2.8V bulbs have silver cases with blue base.















FIGURA 1:1 FIGURE: ACTUAL SIZE	VOLT VOLTS	STRUMENTI INSTRUMENTS	TIPO TYPE	CODICE CODE
	2,8	Oftalmoscopio tascabile - alogena <i>Halogen lamp for Pocket Ophthalmoscope</i>	1011-P-7050	NT1P5282
	3,5	Oftalmoscopio alogeno, Cheratoscopio <i>Halogen Ophthalmoscope, Keratoscope</i>	1011-P-7034	NT1P5055
	2,8	Oftalmoscopio alogeno, Cheratoscopio <i>Halogen Ophthalmoscope, Keratoscope</i>	1011-P-7042	NT1P5223
	6	1,2A per Oftalmoscopio indiretto <i>1,2A for indirect Ophthalmoscope</i>	1012-P-5022	NT2P5022
 scala 1:2	6	Alogena per "FISON" <i>Halogen lamp for "FISON"</i>	1012-P-7001	NT2P5014
	6	Alogena 4,8W per Oftalmoscopio indiretto ad occhiale <i>Halogen lamp 4,8W for spectacle indirect Ophthalmoscope</i>	1012-P-7002	NT2P5102
	6	10W Oftalmoscopio indiretto "ALL PUPILL" e Vantage <i>10W Indirect Ophthalmoscope "ALL PUPILL" and Vantage</i>	1012-P-7003	NT2P5110
	3,6	Oftalmoscopio grandangolo TWIN MAG <i>Wide angle TWIN MAG ophthalmoscope</i>	1012-P-7004	NT2P7004
	6	Oftalmoscopio indiretto VANTAGE PLUS <i>Indirect ophthalmoscope VANTAGE PLUS</i>	1012-P-7006	NT2P7006
	2,8	Retinoscopio tascabile (ad esaurimento) <i>Pocket Retinoscope (selling out)</i>	1013-P-5075	NT3P5075
	2,8	Alogena per Retinoscopio spot (20 mm) (ad esaurimento) <i>Halogen lamp for spot retinoscope (20 mm) (selling out)</i>	1013-P-7000	NT3P5083
	3,5	Alogena per Retinoscopio spot (20 mm) (ad esaurimento) <i>Halogen lamp for spot Retinoscope (20 mm) (selling out)</i>	1013-P-7001	NT3P5091














<b>FIGURA 1:1</b> <i>FIGURE: ACTUAL SIZE</i>	<b>VOLT</b> <i>VOLTS</i>	<b>STRUMENTI</b> <i>INSTRUMENTS</i>	<b>TIPO</b> <i>TYPE</i>	<b>CODICE</b> <i>CODE</i>
	2,8	Alogena per Retinoscopio a striscia (20 mm) (ad esaurimento) <i>Halogen lamp for streak Retinoscope (20 mm) (selling out)</i>	1013-P-7002	<b>NT3P5104</b>
	3,5	Alogena per Retinoscopio a striscia (20 mm) (ad esaurimento) <i>Halogen lamp for streak Retinoscope (20 mm) (selling out)</i>	1013-P-7003	<b>NT3P5112</b>
	2,8	Retinoscopio spot Vista (35 mm) <i>Spot Retinoscope Vista (35 mm)</i>	1013-P-7006	<b>NT3P5155</b>
	3,5	Retinoscopio spot Vista (35 mm) <i>Spot Retinoscope Vista (35 mm)</i>	1013-P-7007	<b>NT3P5163</b>
	2,8	Retinoscopio a striscia Vista (35 mm) <i>Streak Retinoscope Vista (35 mm)</i>	1013-P-7008	<b>NT3P5171</b>
	3,5	Retinoscopio a striscia Vista (35 mm) <i>Streak Retinoscope Vista (35 mm)</i>	1013-P-7009	<b>NT3P5198</b>
	6	Tonometro PULSAIR 2000 e 3000 <i>2000 and 3000 Tonometer PULSAIR</i>	1024-P-7156	<b>NT4P5062</b>
	3,5	Alogena per Ootoscopio medico (ad esaurimento) <i>Halogen lamp for medic Ootoscope (selling out)</i>	1015-P-7058	<b>NT5P5116</b>
	2,8	Alogena per Ootoscopio medico (ad esaurimento) <i>Halogen lamp for medic Ootoscope (selling out)</i>	1015-P-7066	<b>NT5P5175</b>
	2,5	Torcia tascabile <i>Mes lenscapped</i>	1016-P-5003	<b>NT6P5003</b>
	6	Oftalmoscopio indiretto ad occhiale spectra - 6W <i>Indirect ophthalmoscope spectra spectacle - 6W</i>	1299-P-7075	<b>NT99P7075</b>













<b>FIGURA 1:1</b> <i>FIGURE: ACTUAL SIZE</i>	<b>VOLT</b> <i>VOLTS</i>	<b>WATT</b> <i>WATTS</i>	<b>STRUMENTI</b> <i>INSTRUMENTS</i>	<b>TIPO</b> <i>TYPE</i>	<b>CODICE</b> <i>CODE</i>
	3	1,5	Manico per Oftalmoscopio B13, B123, B12345 e Otoscopio <i>Hand Ophthalmoscope model B13, B123, B12345 and Ootoscope</i>	L - 01	NZL01
	6	3	Manico per Oftalmoscopio B12 e Super Magnifying Oftalmoscopio C, Eutiscopio <i>B12 hand Ophthalmoscope, Super Magnifying Ophthalmoscope C, Euthyscope</i>	L - 02	NZL02
	4	3,6	Retinoscopio a Striscia R1 e RX1 <i>Streak Retinoscopio model R1 and RX1</i>	L - 04	NZL04
	3	0,6	Retinoscopio Spot 3R1-SP e RX-3SP <i>Spot Retinoscopio model 3R1 and RX-3SP</i>	L - 05	NZL05
	3	0,9	Otoscopio a-No. 4, laringoscopio a-No. 5, transilluminatore MS - II <i>Otoscope a-No. 4, laryngoscope a-No. 5, transilluminator MS-II</i>	L - 06	NZL06
	12	7,2	Flash Booster, Eutiscopio <i>Flash Booster, Euthyscope</i>	L - 08	NZL08
	6	12	Contacto - Gauge CG <i>Contacto - Gauge CG</i>	L - 11	NZL11
	3	0,9	Medicospot MS <i>Medicospot MS</i>	L - 12	NZL12
	4	0,3	Fissazione luminosa delle lampade a fessura SL-J, SL-CX, CT-R <i>Fixation light of Slit Lamp SL-J, SL-CX, CT-R</i>	L - 13	NZL13
	3	1	Retinoscopio a Striscia 3R1-streak <i>Streak Retinoscope 3R1-streak</i>	L - 15	NZL15
	3,6	1,8	Oftalmoscopio <i>Ophthalmoscope</i>	L - 21	NZL21
	4	2,6	Retinoscopio a Striscia NR, NR-RC, NR-RP, RX-2, RX-RCII, RX-RP <i>Streak Retinoscope NR, NR-RC, NR-RP, RX-2, RX-RCII, RX-RP</i>	L - 27	NZL27
	4	2,6	Retinoscopio a Striscia NR-3, NR-3A, RX-3, RX3A <i>Streak Retinoscope NR-3, NR-3A, RX-3, RX3A</i>	L - 28	NZL28
 scala 1:2	4	2,5	Oftalmoscopio Alogeno BX-12, BX123, BX-RC, BX-RP, Bxa-12-123, Bxa-RCII, Bxa-RP <i>Halogen Ophthalmoscope BX-12, BX123, BX-RC, BX-RP, Bxa-12-123, Bxa-RCII, Bxa-RP</i>	L - 29	NZL29

<b>FIGURA 1:1</b> <i>FIGURE: ACTUAL SIZE</i>	<b>VOLT</b> <i>VOLTS</i>	<b>WATT</b> <i>WATTS</i>	<b>STRUMENTI</b> <i>INSTRUMENTS</i>	<b>TIPO</b> <i>TYPE</i>	<b>CODICE</b> <i>CODE</i>
 scala 1:2	3	2	Oftalmoscopio Alogeno BX-13, BX13A, BX123, BX-12345, BXa-13, BXa-13A, BXa-123, BXa-12345 <i>Halogen Ophthalmoscope BX-13, BX13A, BX123, BX-12345, BXa-13, BXa-13A, BXa-123, BXa-12345</i>	L - 30	<b>NZL30</b>
	3,6	1,6	Retinoscopio Spot Ricaricabile <i>Rechargeable Spot Retinoscope</i>	L - 32	<b>NZL32</b>
	3,6	1,8	Retinoscopio Spot Ricaricabile, Laringoscopio <i>Rechargeable Spot Retinoscope, Laryngoscope</i>	L - 33	<b>NZL33</b>
	3	0,9	Doctor Light DL <i>Doctor Light DL</i>	L - 37	<b>NZL37</b>
	3	2	Oftalmoscopio Tascabile Alogeno GH, a-Jr <i>Halogen Pocket Ophthalmoscope GH, a-Jr</i>	L - 38	<b>NZL38</b>
	6	10	Oftalmoscopio Binoculare Indiretto Alogeno IO-H <i>Halogen Binocular Indirect Ophthalmoscope IO-H</i>	L - 41	<b>NZL41</b>
	6	10	Oftalmoscopio Binoculare Indiretto Alogeno IO-a <i>Halogen Binocular Indirect Ophthalmoscope IO-a</i>	L - 51	<b>NZL51</b>
	10	15	Britescope BS - II <i>BS - II Britescope</i>	L - 52	<b>NZL52</b>
	3	1	Set Otoscopio Jr. B <i>Otoscope Set Jr. B</i>	L - 55	<b>NZL55</b>
	3	2	Otoscopio F.O. FD-AC, AC-R, Illuminatore FD-AD, AD-R <i>Fiber Otoscope F.O. FD-AC, AC-R, Fiber Illuminator FD-AD, AD-R</i>	L - 57	<b>NZL57</b>
	6	5,8	Britescope BS-Jr <i>Britescope BS-Jr</i>	L - 60	<b>NZL60</b>
	3	1,7	Set diagnostico, BXa - 12345FO, Bxa - 134FO, otoscopio fiber FO, a - 34FO <i>Diagnostic set, BXa-12345FO, Bxa - 134FO, fiber otoscope FO, a-34FO</i>	L - 69	<b>NZL69</b>


















<b>FIGURA 1:1</b> <i>FIGURE: ACTUAL SIZE</i>	<b>VOLT</b> <i>VOLTS</i>	<b>STRUMENTI</b> <i>INSTRUMENTS</i>	<b>CODICE</b> <i>CODE</i>
	2,5	Otoscopi, Anoscopi, Transilluminatori <i>Otosopes, Anoscopes, Transilluminators</i>	<b>WA00200</b>
	2,5	Sigmoidoscopio 30800 (ad esaurimento) <i>30800 Sigmoidoscope (selling out)</i>	<b>WA00600</b>
	2,5	Retinoscopi, Oftalmoscopi 10600, 121, 110, 17600, 17610 <i>10600, 121, 110, 17600, 17610 Retinoscopes, Ophthalmoscopes</i>	<b>WA00900</b>
	2,5	Transilluminatore 43100 (fuori produzione - ad esaurimento) <i>43100 Transilluminator (delisted - selling out)</i>	<b>WA01100</b>
	4,65	Alogena - Oftalmoscopio binoculare indiretto <i>Halogen lamp for binocular, indirect Ophthalmoscope</i>	<b>WA01200</b>
	2,5	Abbassalingua, Penne luminose 28000, 76000, 77700 <i>28000, 76000, 77700 tongue blade holder, Pen light</i>	<b>WA01300</b>
	2,5	Penna luminosa 77800 <i>77800 Pen light</i>	<b>WA01400</b>
	2,5	Retinoscopio 17700 (originali WA01800 - fuori produzione) <i>17700 Retinoscope (original WA01800 - delisted)</i>	<b>X-01.88.058</b> <b>(WA01800)</b>
	7	Lampada frontale 45003 <i>45003 Head light</i>	<b>WA02000</b>
	2,5	Oftalmoscopio obsoleto 11500 (ad esaurimento) <i>11500 Ophthalmoscope (selling out)</i>	<b>WA02100</b>
	2,5	Otoscopio 22900 (ad esaurimento) <i>22900 Otoscope (selling out)</i>	<b>WA02200</b>
	2,5	Oftalmoscopio 12900 - obsoleto (ad esaurimento) <i>12900 Ophthalmoscope - delisted (selling out)</i>	<b>WA02300</b>
	6	Lampada frontale 46003 <i>46003 Head light</i>	<b>WA02500</b>
	6	Alogena per Lampada frontale <i>Halogen lamp for Head light</i>	<b>WA02600</b>














<b>FIGURA 1:1</b> <i>FIGURE: ACTUAL SIZE</i>	<b>VOLT</b> <i>VOLTS</i>	<b>STRUMENTI</b> <i>INSTRUMENTS</i>	<b>CODICE</b> <i>CODE</i>
	2,5	Oftalmoscopio tascabile <i>Pocket Ophthalmoscope</i>	<b>WA02800</b>
	2,5	Otoscopio tascabile <i>Pocket Otoloscope</i>	<b>WA02900</b>
	3,5	Alogena per Oftalmoscopi e Retinoscopi 11710, 12400, 47300, 11600, 11605, 11610, 18000 <i>Halogen lamp for 11710, 12400, 47300, 11600, 11605, 11610, 18000 Ophthalmoscopes and Retinoscopes</i>	<b>WA03000</b>
	3,5	Alogena per Otoscopi e Transilluminatori 2000, 25020, 21700, 20200, 20202, 25000... <i>Halogen lamp for 2000, 25020, 21700, 20200, 20202, 25000... Otoloscopes and Transilluminators</i>	<b>WA03100</b>
	2,5	Alogena per Oftalmoscopi 11511, 11500 <i>Halogen lamp for 11511, 11500 Ophthalmoscopes</i>	<b>WA03300</b>
	2,5	Alogena per Otoscopi, Illuminatori, Transilluminatori 21110, 24031, 24011, 24020, 24000, 27200 <i>Halogen lamp for 21110, 24031, 24011, 24020, 24000, 27200 Otoloscopes, Illuminators, Transilluminators</i>	<b>WA03400</b>
	3,5	Alogena per Video-otoscopio <i>Halogen lamp for Video-otoscope</i>	<b>WA03500</b>
	2,5	Dispositivo illuminazione laringoscopia 60500 <i>60500 disposable laryngoscope illuminator</i>	<b>WA03600</b>
	3,5	Alogena per Retinoscopio a fessura 18100 <i>Halogen lamp for 18100 Retinoscope</i>	<b>WA03700</b>
	3,5	Alogena HPX per Panoptic 11810-11820 <i>HPX halogen lamp for Panoptic 11810-11820</i>	<b>WA03800</b>
	2,5	Alogena per Oftalmoscopi 12810, 11110 <i>Halogen lamp for 12810, 11110 Ophthalmoscopes</i>	<b>WA03900</b>
 scala 1:2	21	Fonte di luce LX-150 150W alogena <i>LX-150 light source halogen lamp 150W</i>	<b>WA04050</b>
 scala 1:2	12	35W per lampada LS135 <i>35W for LS 135 lamp</i>	<b>WA04430</b>















<b>FIGURA 1:1</b> <i>FIGURE: ACTUAL SIZE</i>	<b>VOLT</b> <i>VOLTS</i>	<b>STRUMENTI</b> <i>INSTRUMENTS</i>	<b>CODICE</b> <i>CODE</i>
 scala 1:2	12	50W per lampada LS 150 <i>50W for LS150 lamp</i>	<b>WA04450</b>
	2,5	Alogena per Oftalmoscopi 11470, 11475, 11400, 11411 (sostituisce la WA04300) <i>Halogen lamp for 11470, 11475, 11400, 11411 Ophthalmoscopes (WA04300 replacement)</i>	<b>WA04400</b>
	3,5	Alogena per Retinoscopio 18010 <i>Halogen lamp for 18010 Retinoscope</i>	<b>WA04500</b>
	2,5	Alogena per Retinoscopio 17710 <i>Halogen lamp for 17710 Retinoscope</i>	<b>WA04600</b>
	2,5	Laringoscopi per bambini McINTOSH/MILLER 634, 664, 674, 680, 684, 690, 692, 694 <i>634, 664, 674, 680, 684, 690, 692, 694 McINTOSH/MILLER Laryngoscopes for children</i>	<b>WA04700</b>
	2,5	Laringoscopi per adulti McINTOSH/MILLER 634, 640, 664, 674, 680, 684, 690, 692, 694, 40415 <i>634, 640, 664, 674, 680, 684, 690, 692, 694, 40415 McINTOSH/MILLER Laryngoscopes for adults</i>	<b>WA04800</b>
	3,5	Oftalmoscopi 11720, 11730, 11620, 11630 - Alogena 11720, 11730, 11620, 11630 Ophthalmoscopes - Halogen lamp	<b>WA04900</b>
	2,5	Alogena per Laringoscopio a fibre ottiche 60813, 60814, 60715, 60815, 60803, 60804 <i>Halogen lamp for 60813, 60814, 60715, 60815, 60803, 60804 Laryngoscope handles</i>	<b>WA06000</b>
	14,5	Lampada da esame 48435 (sostituisce la 04100) <i>48435 Examination lights (in the place of 04100)</i>	<b>WA06100</b>
	3,5	Alogena per Audioscopi 23300, 23000, 23020, 23040 <i>Halogen lamp for 23300, 23000, 23020, 23040 Audioscopes</i>	<b>WA06200</b>
	12	20W alogena per Lampada da esame 44100, 44110 <i>20W halogen lamp for 44100, 44110 Examination light</i>	<b>WA06300</b>
 scala 1:2	12	Alogena per Lampada da esame 44200 <i>Halogen lamp for 44200 Examination light</i>	<b>WA06400</b>

<b>FIGURA 1:1</b> <i>FIGURE: ACTUAL SIZE</i>	<b>VOLT</b> <i>VOLTS</i>	<b>STRUMENTI</b> <i>INSTRUMENTS</i>	<b>CODICE</b> <i>CODE</i>
	2,5	Retinoscopio 17200 <i>17200 Retinoscope</i>	<b>WA07200</b>
	2,5	Alogena per Penna luminosa 76600 <i>Halogen lamp for 76600 Pen light</i>	<b>WA07600</b>
	6	Alogena per illuminatori a manico per: Speculum Vaginale, Anoscopio, Sigmoidoscopio (sostituisce la WA04000) <i>Halogen lamp for handles of: Vaginal Specula, Anoscope, Sigmoidoscope (WA04000 replacement)</i>	<b>WA07800</b>
	12	Lampada da esame 48100 <i>48100 Examination light</i>	<b>WA08100</b>
	3,5	Retinoscopio a fessura 18200 <i>18200 streak Retinoscope</i>	<b>WA08200</b>
	3,5	Retinoscopio a spot 18300 <i>18300 spot Retinoscope</i>	<b>WA08300</b>
	2,5	Retinoscopio tascabile 18400 (ad esaurimento) <i>18400 pocket Retinoscope (selling out)</i>	<b>WA08400</b>
	3,6	Otoscopio binoculare LUMIVIEW e per sistema di ingrandimento a superficie piatta - Flat Surface Magnifier <i>LUMIVIEW Otoscope and for Flat Surface Magnifier</i>	<b>WA08500</b>
 Non in scala No real size	4,6	Per lampada 78812 <i>For 78812 lamp</i>	<b>WA08800</b>
 Scala 1:2	60	24W - fonte di luce MFI - Solarc - 49502 <i>24W - MFI - Solarc - Light source - 49502</i>	<b>WA09500</b>
 Scala 1:2	-	21W - Solarc Lamp	<b>WA09800</b>



<b>FIGURA 1:1</b> <i>FIGURE: ACTUAL SIZE</i>	<b>MOD.</b>	<b>STRUMENTI</b> <i>INSTRUMENTS</i>	<b>TIPO</b> <i>TYPE</i>	<b>CODICE</b> <i>CODE</i>
		Lampada frontale allo xenon "Vision" <i>"Vision" head light</i>	30857	GI30857
		Lampada frontale per "Clar" e "Storz" <i>Head light for "Clar" and "Storz"</i>	31145	A715/2
		Dermatoscopio <i>Dermatoscope</i>	31186	605
		Dermatoscopio "GIMA 2000" <i>"GIMA 2000" Dermatoscope</i>	31188	138X
		Occhiali "Nistagmo" 4V 0,230A 1W <i>"Nistagmo" 4V 0,230A 1W spectacle</i>	31208	155
		Otoscopio "Ifield" <i>"Ifield" Otoscope</i>	31356	GI31356
		Oftalmoscopio "Ifield" <i>"Ifield" Ophthalmoscope</i>	31357	GI31357
	11M	Otoscopio, Oftalmoscopio, riuniti (ad esaurimento) <i>Otoscope, Ophthalmoscope, united (selling out)</i>	31375	GI31375
	12	Otoscopio 3027N - "Gowllands", Oftalmoscopio, Abbassalingua, riuniti <i>3027N Otoscope "Gowllands", Ophthalmoscope, Tongue depressor</i>	31376	GI31376
	13	Oftalmoscopio "May" <i>"May" Ophthalmoscope</i>	31377	GI31377
	38	Otoscopio Veterinaria (ad esaurimento) <i>Veterinary Otoscope (selling out)</i>	31378	GI31378
	310	Otoscopio Standard, mini Otoscopio e Laringoscopi lame da 1 a 4 <i>Standard Otoscope, Otoscope mini, Laryngoscopes blade 1 to 4</i>	31379	GI31379
	910	Laringoscopio Lama 0 <i>Laryngoscope blade 0</i>	31380	167
	611M	Oftalmoscopio "Mayfield" <i>"Mayfield" Ophthalmoscope</i>	31418	GI31418
	638	Oftalmoscopio "Mayfield" <i>"Mayfield" Ophthalmoscope</i>	31419	GI31419

<b>FIGURA 1:1</b> <i>FIGURE: ACTUAL SIZE</i>	<b>MOD.</b>	<b>STRUMENTI</b> <i>INSTRUMENTS</i>	<b>TIPO</b> <i>TYPE</i>	<b>CODICE</b> <i>CODE</i>
		Oftalmoscopio "Parker" <i>"Parke" Ophthalmoscope</i>	31446	GI31446
		Otoscopio "Parker" 2,5V; alcuni modelli Riester, Kawe, Timesco <i>"Parke" 2,5V Otoloscope, some models Riester, Kawe, Timesco</i>	31449	GI31449
		Otoscopi Gimalux - "Mini" Riester e Kawe <i>Gimalux - "Mini" Riester e Kawe otoscopes</i>	31478	137X
		Alogena per otoscopio e dermatoscopio <i>Halogen lamp for otoscope and dermatoscope</i>	31479	GI31479
		3,5V per otoscopio "RI-Former" <i>"RI-Former" otoscope - 3,5V</i>	31825	GI31825
		3,5V per oftalmoscopio "RI-Former" <i>"RI-Former" ophthalmoscope - 3,5V</i>	31826	GI31826
		Laringoscopio tradizionale "Miller" - per bambini <i>"Miller" Conventional Laryngoscope - for children</i>	34324	167
		Laringoscopio tradizionale "Mc-Int" - per adulti <i>"Mc-Int" Conventional Laryngoscope - for adults</i>	34325	06401
		Laringoscopio fibre ottiche per lame (ad esaurimento) <i>F.O. Laryngoscope for blades (selling out)</i>	34350	GI34350
		Green Spec	34430	0661
		Green Spec Pro (disponibile solo lampadina senza supporto) <i>Green Spec Pro (Only the bulb, without support)</i>	34445 fuori prod. delisted	137X
		Laringscopio F.O. Gima, Kawe, Truphateck <i>Gima, Kawe, Truphateck F.O. laryngoscope</i>	34485	GI34485
		Otoscopio "Fibrex" (ad esaurimento) <i>"Fibrex" Otoloscope (salling out)</i>	fuori prod. delisted	GI43358

<b>FIGURA 1:1</b> <i>FIGURE: ACTUAL SIZE</i>	<b>VOLT</b> <i>VOLTS</i>	<b>STRUMENTI</b> <i>INSTRUMENTS</i>	<b>TIPO</b> <i>TYPE</i>	<b>CODICE</b> <i>CODE</i>
	2,7	Otoscopio UNI I, UNI III, ECONOM, Abbassalingua, Otoscopio SPECULIGHT <i>Otosopes model UNI I, UNI III, ECONOM, Tongue Depressors, SPECULIGHT Otoscope</i>	10421	114
	2,7	Otoscopio UNI II, UNI III, ECONOM <i>Otosopes model UNI I, UNI III, ECONOM</i>	10424	004
	2,7	Otoscopio PEN-SCOPE <i>PEN-SCOPE Otoscope</i>	10472	A648
	2,7	Oftalmoscopio PEN-SCOPE e Varifocal <i>PEN-SCOPE Ophthalmoscope and Varifocal</i>	10473	184
	2,5	Otoscopio, Abbassalingua, Specolo Nasale, Illuminatore Braccio ricurvo RI-SCOPE, RI-STAR, RI-MINI <i>Otoscope, Tongue Depressor, Nasal Speculum, Illuminators model RI-SCOPE, RI-STAR, RI-MINI</i>	10600	600
	3,5	Otoscopio, Abbassalingua, Specolo Nasale, Illuminatore Braccio ricurvo RI-SCOPE <i>Otoscope, Tongue Depressor, Nasal Speculum, RI-SCOPE Illuminator</i>	10601	601
	2,5	Laringoscopio - manico a fibra ottica, Otoscopio operatorio F.O. <i>F.O. Laryngoscopes, Operating Otoscope</i>	10602	904
	3,5	Laringoscopio - manico a fibra ottica, Otoscopio operatorio F.O. <i>F.O. Laryngoscopes, Operating Otoscope</i>	10603	904A
	2,5	Oftalmoscopio RI-SCOPE, RI-STAR, RI-MINI e Dermatoscopio RI-DERMA <i>Ophthalmoscope model RI-SCOPE, RI-STAR, RI-MINI and RI-DERMA Dermatoscope</i>	10605	605
	3,5	Oftalmoscopio RI-SCOPE, RI-STAR, RI-MINI e Dermatoscopio RI-DERMA <i>Ophthalmoscope model RI-SCOPE, RI-STAR, RI-MINI and RI-DERMA Dermatoscope</i>	10606	606
	3,5	Otoscopio RI-MINI, RI-STAR, RI-SCOPE <i>RI-MINI, RI-STAR, RI-SCOPE otoscope</i>	10607	607
	3,5	Oftalmoscopio RI-SCOPE <i>RI-SCOPE ophthalmoscope</i>	10608	608
	2,5	Retinoscopio RI-SCOPE a fessura <i>RI-SCOPE Streak Retinoscope</i>	10615	RI615
	3,5	Retinoscopio RI-SCOPE a fessura <i>RI-SCOPE Streak Retinoscope</i>	10616	RI616

<b>FIGURA 1:1</b> <i>FIGURE: ACTUAL SIZE</i>	<b>VOLT</b> <i>VOLTS</i>	<b>STRUMENTI</b> <i>INSTRUMENTS</i>	<b>TIPO</b> <i>TYPE</i>	<b>CODICE</b> <i>CODE</i>
	2,5	Retinoscopio RI-SCOPE a spot <i>RI-SCOPE Spot Retinoscope</i>	10620	620
	3,5	Retinoscopio RI-SCOPE a spot <i>RI-SCOPE Spot Retinoscope</i>	10621	621
	2,2	Penne Luminose FORTELUX N 0,18A <i>FORTELUX N 0,18A Pen light</i>	11177	A647
	2,5	NORMALUX, CLAR 55, CLAR 100 <i>NORMALUX, CLAR 55, CLAR 100</i>	11300	A625
	6	NOVALUX, CLAR 55, CLAR 100 - senza manico <i>NOVALUX, CLAR 55, CLAR 100 - without handle</i>	11301	RI6106
	6	Lampada frontale RI-CLAR 55mm, RI-CLAR 100mm <i>Head light model RI-CLAR 55mm, RI-CLAR 100mm</i>	11320	320
 scala 1:2	8	RI-FOCUS, REVERSE, RI-SPECIAL, RI-MAGIC	11321	ALO55
	2,5	Laringoscopi (per bambini) - lama 0 - 1 <i>Laryngoscopes (for children) - blade 0 - 1</i>	11378	169
	2,7	Laringoscopi (per adulti) - lama 2 - 3 - 4 <i>Laryngoscopes (for adults) - blade 2 - 3 - 4</i>	11379	171
	2,5	Laringoscopi (per bambini) - lama 0 - 1 <i>Laryngoscopes (for children) blade 0 - 1</i>	11380	167
	2,7	Laringoscopi (per adulti) - lama 2 - 3 - 4 <i>Laryngoscopes (for adults) - blade 2 - 3 - 4</i>	11381	06401
	2,5	Laringoscopio <i>Laryngoscope</i>	11428	428
	3,5	Laringoscopio <i>Laryngoscope</i>	11429	429
	2,5	Toplight handles	11431	431